

I M A G I N E

Imagine a seamless workflow  
from protein to peptide

ProteCol™ HPLC Columns turn imagination into results




# ProteCol™-G C8 HPLC columns solve the problem of separating large proteins, including membrane proteins.

How does ProteCol™ turn imagination to results?

ProteCol™-G C8 columns combine these key features:

- 1000 Å pore size
- Intermediate polarity C8 phase
- Inert glass lined column hardware

Why choose 1000 Å pore size?



100 kDa molecules adsorbed on the wall of a 1000 Å pore




100 kDa spherical molecules in a 300 Å pore



100 kDa rod-shaped molecules in a 300 Å pore (aspect ratio 2:1)

1000 Å pore size silicas enable large irregular shaped proteins to bind to the bonded phase without restricting access to the pore - compared to 300 Å silicas whose pores are easily blocked by large proteins.

ProteCol™-G C8's large pore size, C8 phase plus column hardware turn imagination into results by enabling this impressive recovery of intact membrane proteins.



Base peak chromatogram of mouse liver membrane proteins separated on a ProteCol™-G C8 HQ1003 300 µm x 100 mm column.


Separation of intact membrane proteins derived from mouse liver with ProteCol™-G C8 1000

## ProteCol™-G turns imagination into results by enabling HPLC for membrane proteins

- Continuity of using HPLC for all separation needs simplifies your workflow
- Facilitates the use of MS
- Eliminate SDS-PAGE from your workflow

Can you imagine being able to separate this many proteins on any other HPLC column?

*Micro-preparative separation of proteins extracted from a mouse liver microsomal preparation<sup>#</sup>*




### Examples of integral membrane proteins separated by ProteCol™-G C8

Fraction No.	Description	No. of predicted Transmembrane Helices
8	Transmembrane emp 24-like trafficking protein	2
20	B-cell receptor-associated protein	3
24	Solute carrier family 4 – anion exchange	12

<sup>#</sup> Three consecutive runs overlaid


\* Collected fractions were digested with trypsin and the resulting peptides were separated in a ProteCol™-G C18 G HQ303 and analysed on a Thermo Scientific linear ion trap MS. Resulting mzXML files were searched against the ENSEMBL mouse data base using GPM-Xt software. Proteins were identified with a protein expectation cutoff value of 10<sup>-10</sup>.

Imagine what your proteomics project can achieve with ProteCol™!



# HPLC Column Selection Guide

Whether your proteomics strategy is top down or bottom up - we have the column for you!



## ProteCol™ HPLC Column Specifications and Part numbers

ProteCol™ Column	Stationary Phase	Pore Size (Å)	Particle Size (µm)	Lining	Pore Volume (ml)	Surface Area (m <sup>2</sup> /g)	Carbon Load (%)	Calculated Bonded Phase (µmol/m <sup>2</sup> )	ID Coverage	length (mm)	Part Number
G C8 HQ1003	C8 HQ	1000	3	Glass	0.8 ± 0.1	25 ± 5	0.7	2.33	2.1 mm	150	250170
G C8 HQ1003	C8 HQ	1000	3	Glass	0.8 ± 0.1	25 ± 5	0.7	2.33	2.1 mm	100	250172
G C8 HQ1003	C8 HQ	1000	3	PEEKsil™	0.8 ± 0.1	25 ± 5	0.7	2.33	300 µm	150	250175
G C8 HQ1003	C8 HQ	1000	3	PEEKsil™	0.8 ± 0.1	25 ± 5	0.7	2.33	300 µm	100	250177
G C8 HQ1003	C8 HQ	1000	3	PEEKsil™	0.8 ± 0.1	25 ± 5	0.7	2.33	150 µm	150	250180
G C8 HQ1003	C8 HQ	1000	3	PEEKsil™	0.8 ± 0.1	25 ± 5	0.7	2.33	150 µm	100	250182
G C18 HQ203	C18 HQ	200	3	Glass	1.0 ± 0.1	200 ± 30	12.6	2.63	2.1 mm	150	250150
G C18 HQ203	C18 HQ	200	3	Glass	1.0 ± 0.1	200 ± 30	12.6	2.63	2.1 mm	100	250152
G C18 HQ203	C18 HQ	200	3	PEEKsil™	1.0 ± 0.1	200 ± 30	12.6	2.63	300 µm	150	250155
G C18 HQ203	C18 HQ	200	3	PEEKsil™	1.0 ± 0.1	200 ± 30	12.6	2.63	300 µm	100	250157
G C18 HQ203	C18 HQ	200	3	PEEKsil™	1.0 ± 0.1	200 ± 30	12.6	2.63	150 µm	150	250160
G C18 HQ203	C18 HQ	200	3	PEEKsil™	1.0 ± 0.1	200 ± 30	12.6	2.63	150 µm	100	250162
G C18 HQ303	C18 HQ	300	3	PEEKsil™	1.0 ± 0.1	100 ± 20	7.4	3.08	2.1 mm	150	250135
G C18 HQ303	C18 HQ	300	3	Glass	1.0 ± 0.1	100 ± 20	7.4	3.08	2.1 mm	100	250132
G C18 HQ303	C18 HQ	300	3	Glass	1.0 ± 0.1	100 ± 20	7.4	3.08	2.1 mm	150	250130
G C18 HQ303	C18 HQ	300	3	PEEKsil™	1.0 ± 0.1	100 ± 20	7.4	3.08	300 µm	100	250137
G C18 HQ303	C18 HQ	300	3	PEEKsil™	1.0 ± 0.1	100 ± 20	7.4	3.08	150 µm	150	250140
G C18 HQ303	C18 HQ	300	3	PEEKsil™	1.0 ± 0.1	100 ± 20	7.4	3.08	150 µm	100	250142
<b>Guard and Trap Columns</b>											
G C8 HQ1003	C8 HQ	1000	3	Glass	0.8 ± 0.1	25 ± 5	0.7	2.33	2.1 mm	10	250019
G C8 HQ1003	C8 HQ	1000	3	PEEKsil™	0.8 ± 0.1	25 ± 5	0.7	2.33	300 µm	10	2222051
G C18 HQ203	C18 HQ	200	3	Glass	1.0 ± 0.1	200 ± 30	12.6	2.63	2.1 mm	10	250021
G C18 HQ203	C18 HQ	200	3	PEEKsil™	1.0 ± 0.1	200 ± 30	12.6	2.63	300 µm	10	2222053
G C18 HQ303	C18 HQ	300	3	Glass	1.0 ± 0.1	100 ± 20	7.4	3.08	2.1 mm	10	250029
G C18 HQ303	C18 HQ	300	3	PEEKsil™	1.0 ± 0.1	100 ± 20	7.4	3.08	300 µm	10	2222054

### AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd  
Toll Free: 1800 800 167  
Tel: +61 (0) 3 9837 4200  
Fax: +61 (0) 3 9874 5672  
Email: support@sge.com

### CHINA

SGE Shanghai Representative Office  
Tel: +86 21 6407 9382  
Fax: +86 21 6407 9386  
Email: china@sge.com

### MIDDLE EAST

SGE Gulf  
Tel: +971 6 557 3341  
Fax: +971 6 557 3541  
Email: gulfsupport@sge.com

### EUROPE

SGE Europe Ltd  
European Head Office  
Toll Free: 00800 2790 8999  
Toll Free Fax: 00800 2626 2609  
Tel: +44 1908 568 844  
Fax: +44 1908 566 790  
Tel France: +33 1 69 29 80 90  
Fax France: +33 1 69 29 09 25  
Tel Germany: +49 (0) 6155 / 60746 0  
Fax Germany: +49 (0) 6155 / 60746 50  
Email: europe@sge.com

### INDIA

SGE Laboratory Accessories Pvt Ltd  
Tel: +91 22 24715896  
Fax: +91 22 24716592  
Email: sgeindia@vsnl.com

### UNITED STATES OF AMERICA

SGE Incorporated  
Toll Free: (800) 945 6154  
Tel: +1 512 837 7190  
Fax: +1 512 836 9159  
Email: usa@sge.com

### JAPAN

SGE Japan Inc  
Tel: +81 45 222 2885  
Fax: +81 45 222 2887  
Email: japan@sge.com

 **SGE Analytical Science**  
www.sge.com