

Now sold under the
Thermo Scientific brand



Determination of Trace Anions in Hydrofluoric Acid, Ammonium Fluoride, and a Buffered Oxide Etchant

INTRODUCTION

Concentrated hydrofluoric acid (HF) is used in the semiconductor and electronics industries, both alone and as one component of a buffered oxide etchant (BOE), to remove oxide layers during device production. This HF must be of high purity, especially with respect to anionic impurities that can damage the semiconductor (e.g., chloride and sulfate). In 1997, Semiconductor Equipment and Materials International (SEMI) specified that HF must have chloride and sulfate concentrations $< 200 \mu\text{g/L}$ and nitrate and phosphate concentrations $< 100 \mu\text{g/L}$.¹ Simply diluting the HF to a concentration that will not overload the anion-exchange column does not allow enough sensitivity to determine the contaminating anions.

To address this challenge, Watanabe and Ishzaki used an ion-exclusion (ICE) column to separate the strong acid anions from the fluoride, diverted the strong acid anions to an anion-concentrator column, eluted the anions from the concentrator column, and then determined them by ion chromatography (IC) with suppressed conductivity detection.² The diversion of part of one separation to a concentrator is sometimes called a heart-cutting technique. This basic method later was used by Wu and Chen and also applied to phosphoric acid.³ A few years later, Kaiser, Rohrer, and Watanabe improved the method by using recently introduced concentrator and anion-exchange columns that were better suited for the application, and exploring the factors that impact method accuracy, reproducibility, and ruggedness.⁴ They applied

their method to phosphoric acid, HF, and glycolic acid. The methods for each acid were detailed in Dionex Technical Notes 44, 45, and 46, respectively.⁵⁻⁷

Since the development of this method for determining strong acid anions in weak acids, there have been a number of improvements in IC. Dionex introduced eluent generation to produce high-purity potassium hydroxide eluents, high-capacity anion-exchange columns and concentrators to use with hydroxide eluents, new suppressors that lower detection limits, new IC hardware (ICS-3000 Ion Chromatography system) that facilitates heart-cutting methods, and software improvements that help produce a calibration curve in the samples to be analyzed. Since the publication by Kaiser et al., there have been some reports of using one or more of these improvements.

In 2002, Wang et al. reported using the method of Kaiser et al. but noted that sulfate was difficult to determine because of a high sulfate blank from the IonPac® ICE-AS6 column, and that the reproducibility for phosphate determinations was not good.⁸ The following year, Wilhelm Blödern of Honeywell reported the application of the same method using the IonPac AC10 concentrator and the IonPac AS10 separator column.⁹ He also reported applying the method to a mixture of 25% HF and 20% ammonium fluoride. More recently, Vermeiren modified the method to use an eluent generator to produce ultrapure KOH, an IonPac AC10 concentrator, and an IonPac AS18 separator; he also replaced the

ICE-AS6 column with the ICE-AS1 column.¹⁰ Using an ultrapure hydroxide eluent improves method sensitivity. The IonPac AS18 high-capacity column is designed for use with hydroxide eluents; therefore, it is well suited for this application. The ICE-AS1 column has a much lower sulfate blank than the typical ICE-AS6 column.

This work, an update to the original TN 45, reports an improved method for determining low concentrations of strong acid anions in HF. The method design also allows determination of low concentrations of strong acid anions in an HF/ammonium fluoride mixture (BOE), and ammonium fluoride. The ICE-AS6 column was replaced by the ICE-AS1 column because internal analysis showed that ICE-AS1 column had consistently low levels of sulfate, whereas there was a wide variation in sulfate concentrations from the ICE-AS6 column. An ICS-3000 system with Chromeleon® Chromatography Data System (CDS) software was used to execute the method using the second pump of the dual-pump module to perform ICE separation, and the standard addition feature of the software to produce the calibration curve. The strong acid anions were captured on a 4 mm IonPac AG11-HC guard column that served as the concentrator, and then separated on a 2 mm IonPac AS11-HC column set. After each injection, the ICE-AS1 column was washed with formic acid. This wash had at least two purposes. First, it removed corrosive HF from the column, sample loop, and other tubing. Second, it produced the acidic environment necessary for determining strong acid anions in ammonium fluoride, which is not acidic. The authors also empirically determined that the wash improved reproducibility for determinations of strong acid anions in HF/ammonium fluoride mixtures. The greater sensitivity delivered by using a hydroxide eluent, eluent generation, and newer generation suppressors allowed determinations to be made with 12% HF rather than 24.5% HF, and a 500 µL injection rather than 750 µL. Only one-third of the amount of HF was injected, compared to what was required for the original version of TN 45.

EQUIPMENT

Dionex ICS-3000 system including:

DP Dual Pump

DC Detector/Chromatography module with dual-temperature zone equipped with an injection valve

AM Automation Manager with 10-port valve

EG Eluent Generator module

AS Autosampler

DXP Pump

Chromeleon software Version 6.80 SP6

REAGENTS AND STANDARDS

Deionized water (DI), Type I reagent-grade, 18 MΩ-cm resistivity or better

50% Hydrofluoric acid (HF, Ajax)

98 to 100% Formic acid (CH₃COOH, Merck)

Ammonium hydrogen difluoride (NH₄F.HF, Ajax)
(referred to as ammonium fluoride)

Sodium fluoride (NaF, Fluka)

Sodium chloride (NaCl, Fluka)

Sodium nitrate (NaNO₃, Fluka)

Sodium sulfate (Na₂SO₄, Fluka)

Disodium hydrogen orthophosphate (Na₂HPO₄, Fluka)

400g/L Sodium hydroxide solution (NaOH, Kanto)

PREPARATION OF REAGENTS AND STANDARDS

Eluent Solution

The eluent generator produces the eluent using the EluGen EGC II KOH cartridge and DI water supplied by the pump, with the eluent concentration controlled by the Chromeleon software. Backpressure tubing must be added to achieve 2000 to 2500 psi backpressure that will allow the EG degasser to function properly. See the *ICS-3000 Ion Chromatography System Operator's Manual* (P/N 065031-03) for instructions on adding backpressure.

25% Formic Acid

Add approximately 250 mL DI water to a 500 mL volumetric flask and carefully transfer 125 mL formic acid to the same volumetric flask. Bring to volume with DI water and mix well.

200 mM Sodium Hydroxide

Dilute 20 g of 40% (w/w) sodium hydroxide with degassed DI water to a final weight of 1000 g in an eluent bottle. Avoid the introduction of carbon dioxide from the air into the aliquot of 40% sodium hydroxide and DI water being used to make the eluent.

Table 1. Concentrations of Calibration Standard Stock Solutions and Volumes of 100 mg/L Stock Standard Solutions Used for Standard Additions

Sample	Anion	Concentration (µg/L)			Volume of 100 mg/L Standard Solution (µL) in Total Volume 100 mL		
		1	2	3	1	2	3
		Hydrofluoric acid	Chloride	20	40	80	20
	Nitrate	20	40	80	20	40	80
	Sulfate	20	40	80	20	40	80
	Phosphate	20	40	80	20	40	80
Ammonium fluoride and simulated BOE	Chloride	20	40	80	20	40	80
	Nitrate	200	400	800	200	400	800
	Sulfate	400	800	1600	400	800	1600
	Phosphate	20	40	80	20	40	80

Standard and Sample Solutions

To do the analysis using the method of standard additions, the calibration standard solutions must be prepared in the sample solution, and the concentration of sample in each solution must be the same. To meet these criteria, mix 2× the desired concentration of the working standard solutions and the sample solution in a 1:1 ratio.

100 mg/L Stock Standard Solutions

Prepare 1000 mg/L stock standard solutions of fluoride, chloride, nitrate, sulfate, and phosphate by dissolving 0.221, 0.165, 0.137, 0.148, and 0.149 g of NaF, NaCl, NaNO₃, Na₂SO₄, and Na₂HPO₄, respectively, in separate 100 mL volumetric flasks with DI water. Prepare 100 mg/L stock standard solutions by diluting 10 mL of 1000 mg/L stock standard solutions to 100 mL.

Calibration Standard Stock Solutions

Prepare standards of 2× the desired final concentration in the sample by diluting 100 mg/L anion stock standard solutions using the volumes listed in Table 1.

Sample Stock Solutions

These solutions are prepared at 2× the desired concentration for analysis.

24% Hydrofluoric Acid

Accurately weigh 26 g DI water into a 100 mL polypropylene bottle. Slowly add 24 mL (27.84 g) of 50% HF into the same bottle.

Table 2. Concentrations of Added Standards in the Spiked Sample Solutions (µg/L)

Sample	Anion	Spiked 1	Spiked 2	Spiked 3
Hydrofluoric acid	Chloride	10	20	40
	Nitrate	10	20	40
	Sulfate	10	20	40
	Phosphate	10	20	40
Ammonium fluoride and simulated BOE	Chloride	10	20	40
	Nitrate	100	200	400
	Sulfate	200	400	800
	Phosphate	10	20	40

Simulated BOE Sample Solution

Dissolve 40 g ammonium fluoride in 60 mL DI water. To 30 mL of this solution, carefully add 5 mL of 50% HF and mix. Dilute 20 mL of this solution to 100 mL with DI water and mix.

Ammonium Fluoride Sample Solution

Dilute 25 mL ammonium fluoride solution prepared for BOE preparation with 75 mL DI water and mix.

Calibration Standard and Sample Solutions

The calibration standard solutions were prepared by mixing the calibration standard stock solutions with sample stock solutions 1:1. The sample solutions with known amounts of added standards are referred to as Spiked 1, Spiked 2, and Spiked 3. The added standard concentrations of each spiked sample are shown in Table 2. Stock sample solutions are mixed 1:1 with DI water and referred to as Unspiked.

IONPAC TRAP COLUMN REGENERATION

The Anion Trap Column-High Capacity (ATC-HC) first must be conditioned; the same procedure is used to regenerate the ATC-HC. Monitoring the blank will indicate when regeneration is necessary. Regeneration typically is necessary on a monthly basis, but frequency will depend upon the quality of DI water and usage rate of the instrument. Increased contamination in the water blank indicates that the IonPac ATC-HC needs to be regenerated. The procedure is as follows:

1. Pump 200 mM sodium hydroxide through the ATC-HC at 1.0 mL/min for 50 min
2. Follow with a rinse of DI water at the same flow rate for 20 min.

CONDITIONS

Ion Exclusion

Column: IonPac ICE-AS1, 9 × 250 mm
(P/N 43197)

Eluent: Deionized water

Trap Column: ATC-HC, 9 × 75 mm
(P/N 59640)

Flow Rate: See Table 4

Ion Chromatography

Analytical Column: IonPac AS11-HC, 2 × 250 mm
(P/N 52961)

Guard Column: IonPac AG11-HC, 2 × 50 mm
(P/N 52963)

Concentrator Column: IonPac AG11-HC, 4 × 50 mm
(P/N 53962)

Eluent Source: EGC II KOH (P/N 58900) with
CR-ATC (P/N60477)

Gradient: See Table 4

Flow Rate: 0.38 mL/min

Sample Volume: 500 µL

Column Temperature: 30 °C

Pressure: ~2100 psi

SRS Current: 30 mA

Detection: Suppressed conductivity, Anion
Self-Regenerating Suppressor
(ASRS® 300) (P/N 64555), 2 mm,
external water mode

SYSTEM PREPARATION

To configure the system, refer to Figure 1 and Table 3, which summarize the types and lengths of tubing required for system configuration. The chromatographic hardware is divided into two parts: the ion-exclusion pretreatment portion with an IonPac ICE-AS1 column, and the IC analysis portion with a 2 mm IonPac AS11-HC column set.

To successfully perform this analysis, it is important to use the same type and length of tubing that is listed in Table 3. Changes in tubing length and types will result in a different (and likely incorrect) “cut” from the ICE-AS1 column being delivered to the AS11-HC concentrator column. Tubing colors correspond to the internal diameters listed in Table 3.




Figure 1. System configuration.

Table 3. Details of Tubing Configuration

Connection Points	Tubing Description	Length (cm)	Remark
Port 1=>Port 4 (Injection valve)	Green 0.75 mm (0.033 in)	110	500 µL sample loop
Port 3 (injection valve)=> Port 1 (10-port valve)	Black 0.25 mm (0.010 in)	37	
Port 10=>ICE inlet (10-port valve)	Black 0.25 mm (0.010 in)	4	
ICE outlet=>Port 7 (10-port valve)	Black 0.25 mm (0.010 in)	33	
Port 6=>AG11-HC outlet (10-port valve)	Black 0.25 mm (0.010 in)	5	
AG11-HC inlet=>Port 3 (10-port valve)	Black 0.25 mm (0.010 in)	5	
Port 2=>Port 8 (10-port valve)	Black 0.25 mm (0.010 in)	4	
Port 5=>Analysis (10-port valve)	Red 0.125 mm (0.005 in)	28	Should be as short as possible

Ion Exclusion Pretreatment

1. Prepare the ATC-HC trap column according to directions in the IonPac Trap Column Regeneration section. Caution: Before the ATC-HC trap column is installed in the system, ensure that sodium hydroxide used for storage or cleaning the ATC-HC trap column is completely rinsed away because the ICE-AS1 column is not compatible with hydroxide eluent.
2. Cut a 110 cm portion of the green polyetheretherketone (PEEK™) tubing to make a 500 µL sample loop and install this loop between Ports 1 and 4 of the injection valve.
3. Cut a 37 cm portion of the black PEEK tubing and connect this tubing between Port 3 of the injection valve and Port 1 of the 10-port valve.

- Install the ICE column by connecting a 4 cm piece of the black PEEK tubing between Port 10 of the 10-port valve and the ICE-AS1 column inlet. Use a 32 cm piece of black PEEK tubing between the ICE-AS1 column outlet and Port 7 of the 10-port valve.
- Cut a 5 cm portion of black PEEK tubing and install between Ports 2 and 8 of the 10-port valve.
- Connect the AS injection-port tubing to Port 5 of the injection valve. The injection-port volume must be calibrated before use.
- Cut appropriate lengths of green PEEK tubing for waste lines and connect those pieces of tubing to Port 6 of the injection valve and Port 9 of the 10-port valve.

IC Analysis

- Prepare the ASRS suppressor by following the *QuickStart Instructions* (Dionex Document 031368-01) included with the *Instructions and Troubleshooting Guide for the ASRS*.
- Install the 2 mm IonPac AG11-HC guard and IonPac AS11-HC column set at Port 5 of the 10-port valve using red tubing. To minimize dead volume, use the shortest possible lengths of tubing and ensure that the ends of tubing are cut flush.
- Install the IonPac AG11-HC 4 mm concentrator using two 5 cm pieces of black PEEK tubing. Connect one end of tubing to each end of the IonPac AG11-HC 4 mm column and connect the free ends to Ports 3 and 6 of the 10-port valve in the correct eluent flow direction (inlet connected to Port 6 and the outlet connected to Port 3).
- Put the 10-port valve in position B so that the 4 mm IonPac AG11-HC column is in-line with the 2 mm IonPac AG11-HC and AS11-HC column set. Check for leaks. The expected system backpressure for these three columns at 0.38 mL/min is ~ 2000 psi.




Figure 2. Ion exclusion chromatography of 20 mg/L chloride and 100 mg/L fluoride.

Determining the Fraction Time

To find the fraction (cut) time, the tubing connections were changed so that the conductivity cell was between the ICE column and IonPac AG11-HC concentrator. The DI water flow rate delivered by pump 2 was set to 0.5 mL/min and the 10-port valve was in position A. The mixture of 20 mg/L chloride and 100 mg/L fluoride was injected onto the ICE column. The separation of chloride and fluoride on the ICE column is shown in Figure 2. Based on this experiment, the first 11 min of the ICE separation were established as the collection period for sample analysis. After establishing the fraction time, the conductivity cell was plumbed in the system as shown in Figure 1.

Table 4. Sample Preparation, Analysis, and ICE Column Cleaning Program Description

Program	Time (min)	Flow Rate Pump 2	Injection Valve Position	10-Port Valve Position	Eluent Concentration (mM)	Remark
Sample Preparation	Init	0.5	Load	A	8	AS autosampler starts loading sample to sample loop.
	0	0.5	Inject	A	8	Start fraction collection by flushing sample to concentrator.
	11	Off	Inject	A	8	End fraction collection by switching off pump 2.
Analysis and ICE Column Cleaning	Init	Off	Load	B	8	AS autosampler starts loading 25% formic acid to the sample loop for ICE column cleaning.
	0	0.8	Inject	B	8	Start ICE column cleaning and sample analysis. The ICE column is cleaned by flushing 25% formic acid from the sample loop to the ICE column.
	7	0.8	Inject	B	8	
	15	0.8	Inject	B	8	
	30	0.8	Inject	B	30	
	35	0.5	Inject	B	30	
	40	0.5	Inject	B	8	

SYSTEM OPERATION

Initial experiments established that ICE column cleaning was required after sample injection in order to extend the basic method of the original TN 45 to BOE containing ammonium fluoride, and ammonium fluoride. To accomplish this, two programs were used for each sample injection (Table 4), and two sample lines were written in the sequence for each sample. The sample preparation program (ICE separation) was executed by the first sample line. The subsequent anion analysis and ICE column cleaning were executed by the second sample line.

Sample Preparation Program

For this program, the 10-port valve was in position A. The sample was loaded into the sample loop by the AS sampler and flushed to the ICE column with DI water delivered by pump 2. The strong acid anions excluded by the ICE column were collected on the IonPac AG11-HC concentrator. The sample fraction collection was stopped by switching off the flow of pump 2 at 11 min.

Analysis and ICE Column Cleaning Program

For this program, the 10-port valve was switched to position B. The excluded anions collected on the IonPac AG11-HC concentrator were eluted and separated on an IonPac AS11-HC column set. Formic acid (25%) was loaded to the sample loop by the AS autosampler and flushed to the ICE column by DI water delivered by pump 2.

DISCUSSION OF THE METHOD

This method addresses the challenge of determining trace concentrations of contaminant ions such as chloride, nitrate, sulfate, and phosphate in a matrix composed of a high concentration of fluoride ions. This is accomplished in two steps: ion exclusion chromatography (ICE) to separate the strong acid anions from the large excess of fluoride, followed by an IC analysis of the strong acid anions automatically collected from the ICE separation.

The ion-exclusion mechanism separates ionized species from nonionized or weakly ionized species. This occurs because of repulsion of the strong acid anions by the negatively charged hydration shell on the stationary phase surface called the Donnan membrane. Figure 2 illustrates the application of this ICE mechanism to the separation of 20 mg/L chloride from 100 mg/L fluoride.

To determine the amount of strong acid anions (e.g., chloride and sulfate) in HF, the method of standard additions was used to account for the sample matrix effect. Before analysis, three different concentrations of standards were added to the same volume of sample and the samples were labeled Spiked 1, Spiked 2, and Spiked 3. This yielded three samples with the anion concentrations equal to the amount in the unspiked sample plus the known amount. After analysis of the unspiked and three spiked samples, the peak areas for all samples were plotted against the concentration of the added standard (zero for the unspiked sample). The Y-axis of this plot yielded the anion concentration of the unspiked sample. All described processes were executed




Figure 3. Determination of anions in unspiked 12% HF.

by Chromeleon software. For more information about sequence and QNT file setup for the method of standard additions, go to the Chromeleon software *help* function and view *standard addition*.

RESULTS AND DISCUSSION

Figure 3 shows the chromatography of 12% HF. The concentrations of chloride, sulfate, and other anions are much higher than observed for 24.5% HF used in the original version of TN 45. For example, Figure 3 in this update shows ~ 17 µg/L chloride in 12% HF, whereas the 24.5% HF sample in the original TN 45 had ~ 8 µg/L chloride. When the original work was performed, a customer provided a sample of high-purity HF used for semiconductor applications, which was not available for this update. Therefore, to judge the required sensitivity for the intended application, the method detection limit (MDL) for each anion was estimated using the standard deviation of found anions in seven replicate sample injections and the Student's *t* value for the 99.5% confidence level. Table 5 displays the estimated MDLs in 12% HF. These values are not as low as determined in the original TN 45, but they easily meet the recommended SEMI specifications for 49% HF. For example, the previous work reported an MDL of 0.64 µg/L chloride for




Figure 4. Overlay of chromatograms of HF and three samples of HF spiked with different concentrations of anions.

24.5% HF, while this update reports 1.23 µg/L for 12% HF. Both measures are well under the SEMI specification of 100 µg/L (1.28 and 5.02 µg/L, respectively). Recall that the current method injects only one-third the amount of HF used for the original method. If more sensitivity is needed, the method can be optimized for more concentrated HF and/or a larger injection volume. The percentage HF and injection volume used were chosen to increase method ruggedness.

Before estimating MDLs, the method was calibrated using standard additions with the added standard concentrations reported in Table 2. Figure 4 shows the chromatography of this experiment for 12% HF, and

Table 5. Anion MDLs for Each High Fluoride Sample

Anion	MDL (µg/L)		
	Hydrofluoric Acid (12%)	Ammonium Fluoride (8.33%)	Simulated BOE (10x Dilution)
Chloride	1.23	0.81	0.71
Nitrate	1.01	2.25	3.17
Sulfate	2.49	5.48	1.81
Phosphate	9.86	4.83	1.14

Table 6. Calibration Results Reported by Chromeleon Software

Analyte	# of Points	% R-Squared		
		Hydrofluoric Acid (12%)	Ammonium Fluoride (8.33%)	Simulated BOE (10× Dilution)
Chloride	4	99.52	99.94	100.00
Nitrate	4	99.89	100.00	99.99
Sulfate	4	99.98	99.97	100.00
Phosphate	4	99.90	99.94	99.85




Figure 5. Determination of anions in ammonium fluoride.




Figure 6. Chromatogram of an unspiked BOE sample (10× dilution).

Table 6 reports the calibration results for all three samples. Calibration is linear for all samples, suggesting that the fraction of the ICE column sent to the concentrator column is not overloading the concentrator column and causing analyte loss. The good peak shapes in chromatograms for the unspiked 8.33% ammonium fluoride (Figure 5) and the unspiked simulated BOE sample (Figure 6) also suggest that overloading is not occurring. As an accuracy check, the authors prepared fresh samples, spiked them with 10 µg/L each of chloride, nitrate, sulfate, and phosphate, and then measured the recoveries.

Table 7. Recovery Results

Sample		Anion			
		Chloride	Nitrate	Sulfate	Phosphate
Hydrofluoric Acid (12%)	Spiked Concentration (µg/L)	10	10	10	10
	Concentration in Sample (µg/L)	17.19	1.52	26.11	29.28
	Concentration in Spiked Sample (µg/L)	26.30	10.51	35.92	40.22
	Recovery (%)	96.73	91.20	99.46	102.4
Ammonium Fluoride (8.33%)	Spiked Concentration (µg/L)	10	100	200	10
	Concentration in Sample (µg/L)	80.89	341.6	1021	3.62
	Concentration in Spiked Sample (µg/L)	92.58	438.4	1224	12.97
	Recovery (%)	101.9	99.28	100.2	95.20
Simulated BOE (10x Dilution)	Spiked Concentration (µg/L)	10	100	200	10
	Concentration in Sample (µg/L)	58.79	227.9	694.9	3.10
	Concentration in Spiked Sample (µg/L)	68.86	324.6	893.2	11.88
	Recovery (%)	100.1	99.01	99.81	90.63

Table 8. Retention Time Reproducibility for 51 Sample Injections

Retention Time	Anion			
	Chloride	Nitrate	Sulfate	Phosphate
Average Retention Time (min)	5.51	10.30	14.09	18.27
RSD (%)	0.42	0.62	0.27	0.07

All four anions showed recoveries greater than 90% in all three samples (Table 7). Phosphate recovery is a good measure of method success. Phosphoric acid is the weakest acid within the group that includes hydrochloric, nitric, and sulfuric acids. Therefore, phosphoric acid will be the first anion affected if the properties of the ICE column change, thereby altering the fraction concentrated from the ICE column so that not all phosphate is captured. It is also possible that a change in the ICE column and fraction concentrated can deliver too much fluoride to the concentrator column, thereby overloading it. Because phosphate is trivalent, its retention and recovery will be impacted before that of the other anions.

To assess method reproducibility, retention time precision was calculated using an experiment that evaluated all three samples. For each of the three sample types, there were 17 injections that included samples spiked with known quantities of anions. Table 8 shows that retention time was reproducible for these 51 injections; no trend toward longer or shorter time was observed.

Precautions

Exercise extreme caution when handling concentrated HF, which can corrode glass containers. Therefore, plastic containers must be used for handling HF. High-density polypropylene, high-density polyethylene, or Teflon containers are recommended and Teflon is preferred. To minimize contaminating ions, containers should be soaked with 17.8 MΩ-cm DI water (or better) for at least 24 h before use. All sample preparations with HF should be done in a fume hood. Work behind a shield and wear a face shield, goggles, and gloves designed for handling HF. Calcium gluconate gel should be available in case of HF contact with exposed skin. Consult your safety officer before working with HF.

The IonPac ICE-AS1 column is not compatible with hydroxide eluent, so sodium hydroxide must be flushed out of the ATC-HC trap column after regeneration or initial startup of the trap column. The backpressure of the ICE-AS1 column should not exceed 1000 psi, as noted in the *IonPac ICE-AS1 Product Manual*. Therefore, set the maximum pressure limit to 1000 psi to avoid column damage.

The success of this application depends on a consistent flow rate of DI water delivered by the pump to the ICE-AS1 column in order to produce the correct fraction on the IonPac AG11-HC concentrator. Periodically recalibrate the pump flow rate to ensure delivery of the appropriate fraction.

Long-term exposure to concentrated HF will damage the ICE-AS1 column and PEEK tubing. Rinsing the ICE-AS1 column and sample loop with 25% formic acid and using only 12% HF should improve column and tubing lifetime, compared to earlier methods; however, the ICE-AS1 column may require replacement after approximately 500 injections. Using more concentrated HF solutions will shorten column lifetime. For best results, replace the sample loop and all tubing that comes in contact with HF at least once a year.

REFERENCES

1. *SEMI International Standards: Semiconductor Equipment and Materials International*, Mountain View, CA, Chemical/Reagents Volume, 1997.
2. Watanabe, K.; Ishizaki, Z. Poster Presentation at the International Ion Chromatography Symposium, Dallas, TX, 1995.
3. Wu, M.; Chen, J. *Micro* **1997**, 15 (1), 31–37.
4. Kaiser, E.; Rohrer, J.; Watanabe, K. *J. Chromatogr., A* **1999**, 850, 167–176.
5. Dionex Corporation, *The Determination of Trace Anions in Concentrated Phosphoric Acid*. Technical Note 44, LPN 1084-01, 2002, Sunnyvale, CA.
6. Dionex Corporation, *The Determination of Trace Anions in Concentrated Hydrofluoric Acid*. Technical Note 45, LPN 1106-01, 2002, Sunnyvale, CA.
7. Dionex Corporation, *The Determination of Trace Anions in Concentrated Glycolic Acid*. Technical Note 46, LPN 1057-01, 2002, Sunnyvale, CA.
8. Wang, K.; Lei, Y.; Eitel, M.; Tan, S. *J. Chromatogr. A* **2002**, 956, 109–120.
9. Blödorn, W. Presentation at Semicon Europa, 2003.
10. Vermeiren, K. *J. Chromatogr. A* **2005**, 1085, 60–65.

ASRS, Chromeleon, and IonPac are registered trademarks of Dionex Corporation.
PEEK is a registered trademark of Victrex, Ltd., Company.

Passion. Power. Productivity.



Dionex Corporation

1228 Titan Way
P.O. Box 3603
Sunnyvale, CA
94088-3603
(408) 737-0700

North America

U.S./Canada (847) 295-7500

South America

Brazil (55) 11 3731 5140

Europe

Austria (43) 1 616 51 25 Benelux (31) 20 683 9768 (32) 3 353 4294
Denmark (45) 36 36 90 90 France (33) 1 39 30 01 10 Germany (49) 6126 991 0
Ireland (353) 1 644 0064 Italy (39) 02 51 62 1267 Sweden (46) 8 473 3380
Switzerland (41) 62 205 9966 United Kingdom (44) 1276 691722

Asia Pacific

Australia (61) 2 9420 5233 China (852) 2428 3282 India (91) 22 2764 2735
Japan (81) 6 6885 1213 Korea (82) 2 2653 2580 Singapore (65) 6289 1190
Taiwan (886) 2 8751 6655

www.dionex.com



LPN 1106-02 PDF 08/10
©2010 Dionex Corporation