

Mixed-Mode, Weak Anion-Exchange, Solid-Phase Extraction Method for the Extraction of Niflumic Acid from Human Plasma

Ken Meadows, Thermo Fisher Scientific, Runcorn, UK

Key Words

SOLA WAX, Accucore, mixed mode, weak anion exchange, reproducibility, matrix clean up, SPE

Abstract

An HPLC-MS/MS method has been developed for the determination of niflumic acid in human plasma. This application note demonstrates the use of Thermo Scientific™ SOLA™ WAX SPE products, which allow mixed-mode extractions to be carried out with strong acids, giving improved removal of matrix interferences and reproducible results. The use of a Thermo Scientific™ Accucore™ HPLC column provided fast and efficient separation without the need for an ultra high pressure system. MS/MS detection was performed on a Thermo Scientific™ TSQ Vantage™ mass spectrometer.

Introduction

Niflumic acid is a drug used for treatment of joint and muscular pain. It is categorized as an inhibitor of cyclooxygenase-2. In experimental biology, it has been employed to inhibit chloride channels. It has also been reported to act on GABA and NMDA channels and to block T-type calcium channels.

Niflumic acid is a strong acid that remains ionized across the whole pH range. As a result, mixed-mode, strong anion exchange solid-phase extraction (SPE), which is often the first choice for clean reproducible extraction of acidic compounds, can be challenging as it is not possible to disrupt the ion exchange mechanism by pH adjustment. Extractions can be achieved by using salts as a competitive counter ion; however, the high concentration of salts in the elution solvent can be undesirable particularly for LC-MS/MS analysis. This limits options for extraction to reversed-phase SPE, which can be difficult due to the polar nature and therefore weak retention of many strong acids. Alternative approaches to SPE such as protein precipitation can achieve high recovery of the compounds of interest. However, these approaches tend to exhibit poor reproducibility and significant modification of ionization compared to SPE as they do little to remove endogenous interferences and minimize matrix effects.

SOLA WAX is a mixed-mode, polymeric, weak anion-exchange product that introduces additional selectivity



into the SOLA SPE range. SOLA WAX polymeric SPE utilizes a hydrophobic backbone functionalized with a primary amine with a pKa of approximately 6. This allows users to control the ionization of the stationary phase by modifying the pH. As a result, it is possible to develop simple and effective SPE methods for strong acids that take advantage of the additional matrix cleanliness mixed-mode extractions can offer and the reproducibility benefits of SOLA SPE. SOLA solid phase extraction products introduce next generation, innovative technological advancements, giving unparalleled performance characteristics compared to conventional SPE, phospholipid, and protein precipitation products.

These include:

- Higher levels of reproducibility
- Reduced sample and solvent requirements
- Higher levels of extract cleanliness
- Increased sensitivity

SOLA SPE plates or cartridges provide significant advantages when analyzing compounds in complex matrices, particularly in high-throughput bioanalytical and clinical research laboratories where reduced failure rates, higher analysis speed, and lower sample/solvent requirements are critical. SOLA products' superior performance gives higher confidence in analytical results and lowers cost without compromising ease-of-use or requiring complex method development.

Accucore HPLC columns use Core Enhanced Technology™ to facilitate fast and highly efficient separations. The 2.6 µm diameter particles are not totally porous, but instead have a solid core and a porous outer layer. The optimized phase bonding creates a series of high coverage, robust phases. Accucore RP-MS uses an optimized alkyl chain length for more effective coverage of the silica surface. This coverage results in a significant reduction in secondary interactions and yields highly efficient peaks with very low tailing. The tightly controlled 2.6 µm diameter of Accucore particles results in much lower backpressures than typically seen with sub-2 µm materials.

Experimental Details

| Consumables | Part Number |
|--|-------------|
| Fisher Scientific™ LC/MS grade water | W/011217 |
| Fisher Scientific LC/MS grade methanol | M/4062/17 |
| Fisher Scientific analytical grade formic acid | F/1900/PB08 |

| Sample Handling Equipment | Part Number |
|---|--------------|
| Thermo Scientific™ HyperSep™ glass block manifold 24 port | 60104-233 |
| Stopcocks for 24 port manifold | 60104-244 |
| Vacuum pump | 60104-241 |
| Thermo Scientific™ Chromacol™ WebSeal™ UltraVap™ High Speed Sample Concentrator | CLS-229070 |
| Thermo Scientific™ eVol™ Sample Dispensing System | 66002-024 |
| Thermo Scientific™ Finnpiette™ F2 pipettor kit | PMP-020-220F |

| Vials and Closures | Part Number |
|---|---------------|
| Wide Open Short Screw Thermo Scientific™ Chromacol™ SureStop™ Vial, Clear with ID Patch, GOLD Grade and Red PTFE/White Silicone Septum, Convenience kit | 2-SVWGKST-CPK |

Sample Pretreatment

100 µL of human plasma diluted 1:1 with 2% phosphoric acid

| Sample Preparation | Part Number | |
|---------------------|---|-----------|
| Compound(s): | Niflumic Acid, niflumic acid D3 (IS) | |
| Matrix: | Human plasma | |
| SPE: | SOLA WAX 10 mg 1 mL cartridge | 60109-005 |
| Conditioning stage: | 500 µL methanol | |
| | 500 µL water | |
| Application stage: | Load sample at 0.5 mL/min | |
| Washing stage: | 500 µL 25 mM ammonium acetate (NH ₄ +CH ₃ COO-) buffer in water | |
| | 500 µL methanol | |
| Elution stage: | 500 µL methanol with 2% ammonium hydroxide (NH ₄ OH) at 0.5 mL/min | |
| | Dry under nitrogen and reconstitute in 100 µL 2% formic acid in water | |

| Separation Conditions | | Part Number |
|---------------------------|---|--------------|
| Instrumentation: | Thermo Scientific™ Dionex™ UltiMate™ 3000 LC system | |
| Column: | Accucore RP-MS 2.6 µm, 50 × 2.1 mm | 17626-052130 |
| Guard column: | Thermo Scientific™ Accucore™ RP-MS Defender™ | 17626-012105 |
| | Thermo Scientific™ Uniguard™ drop-in guard holder | 852-00 |
| Flow rate: | 750 µL/min | |
| Run time: | 3 min | |
| Column temperature: | 30 °C | |
| Injection details: | 2 µL full loop injection | |
| Injection wash solvent 1: | Water | |
| Injection wash solvent 2: | IPA / acetonitrile / acetone (45:45:10, v/v/v) | |
| Mobile phase A: | Water + 0.1% formic acid | |
| Mobile phase B: | Methanol + 0.1% formic acid | |
| Gradient conditions: | Table 1 | |

| Time (min) | %A | %B |
|------------|----|-----|
| 0.0 | 60 | 40 |
| 2.0 | 0 | 100 |
| 2.1 | 60 | 40 |
| 3.0 | 60 | 4 |

Table 1: Gradient conditions

MS Conditions

| | |
|------------------------|-------------------------------|
| Instrumentation: | TSQ Vantage Mass Spectrometer |
| Ionization conditions: | HESI |
| Polarity: | Positive |
| Spray voltage: | 3000 V |
| Vaporiser temperature: | 475 °C |
| Sheath gas pressure: | 50 arb |
| Aux gas pressure: | 60 arb |
| Capillary temp: | 300 °C |
| Collision pressure: | 1.5 mTorr |
| Scan time: | 0.02 s |
| Q1: | 0.7 FWHM |
| Q3: | 0.7 FWHM |

Compound transition details are provided in Table 2.

| Compound | Parent (m/z) | S-Lens (V) | Product (m/z) | Collision Energy (V) |
|-----------------------|--------------|------------|---------------|----------------------|
| Niflumic acid | 283 | 78 | 265.0 | 21 |
| Niflumic acid D3 (IS) | 288 | 90 | 270.1 | 22 |

Table 2: Compound transition details

Data Processing

| | |
|-----------|---|
| Software: | Thermo Scientific™ LCQUAN™ software version 2.6 |
|-----------|---|

Results

Niflumic acid standards, extracted from human plasma, gave a linear calibration curve over the dynamic range of 1 to 1000 ng/mL with an R^2 coefficient of 0.999 (Figure 1 and Table 3). The chromatography of the LCQ sample at 3 ng/mL is shown in Figures 2 and 3. QC samples were analyzed at concentrations of 3, 15, 600, and 800 ng/mL (Table 3). Overspikes (post-extraction fortified blank samples) were analyzed at a concentration of 3 ng/mL and used to calculate the percentage recovery level of niflumic acid at 86.6 % (Table 4).

Replicate extractions ($n=20$) were also made at the 3 ng/mL level to assess the reproducibility of the assay. The % CV was calculated at 6.3% for niflumic acid with no internal standard correction and 4.3% when calculated using the response ratio of the internal standard (Table 5). This data shows that the extraction is reproducible even without the addition of an internal standard. The reproducibility of the assay is significantly better than values quoted in literature for a protein precipitation extraction for a 10 ng/mL sample (13.5% $n=6$).[1]

Matrix effects for niflumic acid were calculated at 7.96% (Table 6) compared to 94% quoted for a protein precipitation method.[1] This highlights the improved interference removal and reduction of matrix effects achieved with a mixed-mode SPE extraction.




Figure 2. Example chromatogram of 3 ng/mL niflumic acid




Figure 3: Niflumic acid linearity over the dynamic range 1–1000 ng/mL

| Standard | Specified Conc | Calculated Conc | % Diff |
|----------|----------------|-----------------|--------|
| S1 | 1.0 | 0.940 | -6 |
| S2 | 10.0 | 9.80 | -2 |
| S3 | 25.0 | 27.0 | 8 |
| S5 | 100.0 | 102 | 2 |
| S6 | 500.0 | 491 | -2 |
| S7 | 750.0 | 724 | -3 |
| S8 | 1000.0 | 1030 | 3 |
| QC L | 15.0 | 16.2 | 8 |
| QC M | 600.0 | 609 | 2 |
| QC H | 800.0 | 826 | 3 |

Table 3: Accuracy data for the calibration range 1–1000 ng/mL

| Standard | Response | Absolute Recovery |
|-------------------------|----------|-------------------|
| Average area response | 199024 | 86.6% |
| Overspike area response | 229655 | |

Table 4: Recovery for niflumic acid at 3 ng/mL

| Standard | Response | Matrix Effect |
|-------------------------|----------|---------------|
| Standard area response | 212708 | 7.9% |
| Overspike area response | 229655 | |

Table 6. Matrix effects at 3 ng/mL

| Compound | % CV (n=20) |
|------------------|-------------|
| Niflumic acid | 6.3% |
| Niflumic acid D3 | 7.4% |
| Response ratio | 4.3% |

Table 5: Precision data at 3 ng/mL (n=20)

Conclusion

- SOLA WAX chemistry allows for the fast and easy extraction and quantification of niflumic acid from human plasma.
- Extraction recovery was 86.6%.
- SOLA WAX gave excellent precision for the extraction with % CV (n=20) less than 6.3% even without internal standard correction.
- SOLA WAX achieved low matrix suppression effects at less than 7.9%.
- Accucore RP-MS columns achieved a fast run time of less than 3 minutes.

Reference

- [1] Kang, W., Determination of Talniflumate and Niflumic Acid in Human Plasma by Liquid Chromatography–Tandem Mass Spectrometry, *Analytical Science*, 2009, 2, 571 -574.

thermoscientific.com/chromatography

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

USA and Canada +1 800 332 3331
France +33 (0)1 60 92 48 34
Germany +49 (0) 2423 9431 20 or 21
United Kingdom +44 (0)1928 534110
Japan +81 3 5826 1615

China +86 21 68654588 +86 10 84193588
+86 20 83145199 800 810 5118
India +91 22 6742 9494 +91 27 1766 2352
Australia 1 300 735 292 (free call domestic)
New Zealand 0800 933 966 (free call domestic)
All Other Enquiries +44 (0) 1928 534 050

Technical Support
North America +1 800 332 3331
Outside North America +44 (0) 1928 534 440

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific