

[The Importance of a Perfect Connection in Capillary HPLC

Richard Henry, Hans-Jürgen Wirth, Jana Habsuda,
Peter Dawes, Ern Dawes
SGE International Pty Ltd, 7 Argent Place,
Ringwood, Victoria 3134, Australia

© Copyright 2004 SGE International Pty. Ltd. This presentation cannot be copied or altered in part or full without the expressed written permission of SGE International



Abstract

In small-bore HPLC, especially in capillary column dimensions, the quality of the separation can be greatly reduced by factors outside of the analytical column. While the influence of the end-user on fundamental pump and detector design may be limited, the choice of connection tubing and the way connections are made by the chromatographer can have a very significant impact on data quality.

The objective of capillary LC is the quantification of very small sample amounts, such as a low abundance impurity in a pharmaceutical product or a trace amount of protein or peptide in a biological sample.

Both adsorption of sample components on tubing walls or other materials and excessive band-spreading or mixing in open flow channels can lead to significant problems in capillary dimensions. This paper will focus on the latter issue.

[Small-Bore LC Column Evolution

- It has been 100 years since Tswett published his pioneering LC work on plant pigment separation by what he called “column chromatography”.
- Main small-bore column advantages of greater mass sensitivity and easier interface to certain detectors such as mass spectrometry were documented in the early 1980s by Scott (1) and others.
- Practical applications of small-bore columns, especially below 2mm ID, have been rather slow to develop due to the very limited availability of optimized LC instruments, columns and accessories.


Some Operating Characteristics of Small-Bore Columns

ID (mm)	Typical Flow Rate (mL/min)	Typical Flow Rate (uL/min)	Mass Sensitivity	Solvent Consumption
4.6	1.000000	1000	1	1
4	0.756144	760	1.3	0.76
2	0.189036	200	5.3	0.2
1	0.047259	47	21	0.047
0.5	0.011815	12	85	0.012
0.3	0.004253	4.2	235	0.0042
0.15	0.001063	1.1	940	0.0011
0.075	0.000266	0.27	3762	0.00027


It is a misconception that small-bore LC columns yield greater efficiency than larger ID columns of the same length; more plates can only be realized when small-bore columns of greater length can be made and other components can be optimized.

Brewed Coffee

ProteCol C18, 3 μ m, 300A



PAHs on 300 μ m ID ProteCol™ LC Capillary Column



Column: 300 μ m x 150mm
ProteCol C18, 3 μ m, 120Å
80% AcCN : 20% water
4.0 μ l/min


16 polycyclic aromatic hydrocarbons (EPA 610)

50:1 flow split after injector.

Peak volumes defined by 4s can be less than 1 μ L in capillary LC columns; therefore, volumes of other components must be much smaller than 1 μ L


Peak or Band Dispersion

- Peaks typically show Gaussian or random distribution defined by the standard deviation, s , expressed in time, volume or distance units.
- Dilution or dispersion by mobile phase occurs as solutes pass through the chromatographic system.
- Dispersion can occur both inside and outside of the column bed. Separation gained in the column can be lost in connectors, especially after the column.



Gaussian Peak Shape

$$\sigma_{\text{col}} = V_0 (1 + k) / N^{1/2} \text{ (dispersion in a packed column bed)}$$



Peak volume is often estimated as $4s$ at 13.4 % of peak height. It is proportional to column bed volume, V_0 .

Peak volume is determined by the effects of the column bed (equation) **plus the effects of other system components.**

Broad distribution destroys resolution and sensitivity.

Dispersion in Column Bed^{2,3}

$$\sigma^2 = V_o^2(1 + k')^2/N$$

where variance has units of μL^2 when

V_o = mobile phase volume of column in μL
(unretained peak retention volume)


k' = capacity factor of the peak

N = number of theoretical plates of column

Small geometry, short retention and high efficiency favor low dispersion in columns.

Dispersion in Open Connectors^{2,3}

$$\sigma^2 = 1.36 \times 10^{-3} d_t^4 L_t F/D$$



Velocity at the wall is essentially zero under laminar flow conditions. Small inside diameter, short length, low flow and fast diffusion favor low dispersion in connection tubes and accessories. Larger molecules show greater dispersion ($1/D$).


Total Peak (Band) Dispersion Expressed as Variance Equation^{3,4}

$$\sigma_{\text{tot}}^2 = \sigma_{\text{col}}^2 + \underbrace{\sigma_{\text{inj}}^2 + \sigma_{\text{det}}^2 + \sigma_{\text{conn}}^2}_{\text{System Constant } (\sigma_{\text{sys}})}$$

- Component design determines the System Constant; they often can be optimized by the user for different column situations.
- Sometimes an additional variance term is added to account for thermal mismatch of solvent temperature as sample passes from ambient into a heated or cooled column bed.


Example of a Capillary LC Flow Path

Small-volume components are critical to success




Effect of Tubing Dispersion on 100x2mm Column


.010 inch I.D.



5-cm long






20-cm long




80-cm long

.005 inch I.D.




Effect of Tubing Location on Dispersion³

isocratic



gradient



Tubing after the column always disperses the sample

Tube at detector end

Tube at injector end

Gradient Elution Advantages


- Isocratic elution creates increasing dispersion as a function of efficiency (N), void volume (V_0) and retention (k), causing lower sensitivity for more retained solutes.
- Gradient elution focuses sample components at the column inlet and creates uniform dispersion for all solutes by reducing or eliminating the retention factor aspect.
- In the past, gradients have not been practical for trace analysis due to solvent-related detector noise and drift.
- MS with volatile mobile phases does not have this solvent limitation, so gradient experiments can be conducted at high sensitivity if reliable solvent programs can be generated at low Capillary LC flow rates.

Optimizing Capillary LC Components

The ProteCol™ System

- Integrated approach to Capillary LC incorporating columns, unions, tubing, splitters, filters and accessories.
- All components perfectly matched for minimum dispersion volume.
- **Featuring PEEKsil™ fused silica lined PEEK tubing.**
- Unique, convenient designs virtually eliminate the dispersive effects of components having drilled holes.
- Create a complete column system or improve a system with ProteCol components.

Fused Silica Shows Advantages



One of four drawing towers located at SGE Melbourne, AU

Fused Silica Lined Tubing Benefits⁵



PEEK Tubing:

- Often produces a rough-cut
- Potential dead volume
- Structural integrity greatly affected by solvent, temperature and pressure



PEEKsil Tubing:


- Square cut and polished ends
- No dead volume
- Silica liner protects structural integrity of polymer sheath

A microscope or good magnifying glass is an important tool for inspecting connectors in Capillary LC.

Volume Contribution of Standard PEEKsil™ Connection Tubing Lengths

Internal Tubing Volume in μL							
		Tubing length [cm]					
		5	10	15	25	50	
Tubing ID [μm]	25	0.00	0.00	0.01	0.01	0.02	← 0.001" ID
	50	0.01	0.02	0.03	0.05	0.10	← 0.002" ID
	75	0.02	0.04	0.07	0.11	0.22	← 0.003" ID
	100	0.04	0.08	0.12	0.20	0.39	← 0.005" ID
	150	0.09	0.18	0.27	0.44	0.88	
	175	0.12	0.24	0.36	0.60	1.20	
	200	0.16	0.31	0.47	0.79	1.57	
	300	0.35	0.71	1.06	1.77	3.53	← 0.010" ID





ProteCol™ Components



Columns, guards,
traps, unions, filters
and tubing.


[ProteCol™ Components Employ a Unique Ferrule Concept

- Keeps system dead volume and total volume to an absolute minimum.
- Allows tubing or column connections without adding volume from the fitting body.



ProteCol™ Unions

- Butt tubing connections avoid dead volume and extra volume.
- Each ProteCol Union has a different ferrule.
- Flats for wrench tightening at higher pressures.
- Knurls allow finger tightening at lower pressures.




Industry-wide dimension standards are not yet in place for tubing and fittings. SGE tubing, fittings and ferrules are closely controlled for tight fit and reliable high pressure use.

Influence of Multiple 50 μm ID Connections on Peak Shape



ProteCol C18 100mm x 300 μm Column extended with a total of 40cm PEEKsil tubing (50 μm ID) incorporating 1 to 8 zero volume connections.

Flow rate: 4.0 $\mu\text{l}/\text{min}$
Mobile phase: 60% ACN in water
Naphthalene peak shown




ProteCol™ Capillary Analytical Column

- No endfittings required.
- Efficient butt connections
- Rigid 1/16" and flexible 1/32" versions available.




[ProteCol™ Column Endfitting⁶

Porous element




ProteCol™ Capillary Guard Column



- Bed is 1 cm length.

Test Chromatograms on ProteCol™ Capillary Columns



Column: 300 μ m x 150mm
ProteCol C18, 3 μ m, 120Å
60% AcCN : 40% water
4.0 μ l/min

1 Pyridine
2 Methyl Benzoate
3 Phenetole
4 Naphthalene



50:1 flow split between
injector and column.

Current Phases for the ProteCol™ Column Range

Phase	dp	ligand	pore size	carbon load	surface area	pore volume
ProteCol C18-120-3	3 µm	C18 endcapped	120 Å	17%	350m ² /g	1.1ml/g
ProteCol C18-120-5	5 µm	C18 endcapped	120 Å	17%	350m ² /g	1.1ml/g
ProteCol C18-300-3	3 µm	C18 endcapped	300 Å	9%	100m ² /g	0.9ml/g
ProteCol C18-300-5	5 µm	C18 endcapped	300 Å	9%	100m ² /g	0.9ml/g
ProteCol C4-300-3	3 µm	C4	300 Å	3%	100m ² /g	0.9ml/g
ProteCol NH2-300-3	3 µm	-NH2	300 Å	2%	100m ² /g	0.9ml/g

- ProteCol C18-120 phases are optimized for separation of small molecules and are highly deactivated.
- ProteCol C18-300 phases are optimized for separation of biomolecules such as peptides and proteins.

Van Deemter Plots for ProteCol™ LC Capillary Columns



- ProteCol C18-120-3 (3µm; 120Å)

ProteCol™ Capillary Trap Column



- Allows preconcentration of a sample before injection onto the analytical column.
- Short, flexible design allows easy, efficient valve installation.
- 10mm length standard

ProteCol™ Capillary LC Filtering Connector




1/32" OD



1/16" OD

- Two filter versions protect Capillary LC columns (2um and 50 micron tubing is standard)


Affect of Filtering Connector on Column Performance



Column: 300 μ m x 150mm
ProteCol C18, 3 μ m, 120Å
60% AcCN : 40% water
4.2 μ l/min

Test Mix:
Pyridine
Methyl Benzoate
Phenetole
Naphthalene

Back Pressure Contribution of the Capillary Filtering Connector



The pressure drop is entirely across the tubing, not the filter element. A larger ID Micro version of the filtering connector is available for 1 and 2mm ID columns.

The 2 μ m porosity filter connected to 50 μ m ID tubing of 150mm total length gave a slope of only 2 psi per 1 μ L/min flow, which should be ideal for columns of 1mm ID or smaller.

Backpressure in Capillaries

The backpressure generated in a capillary is described by the Poiseuille's equation:

$$\Delta p = \frac{F \cdot 8\eta \cdot l}{\pi \cdot r^4}$$

Δp = pressure drop


η = viscosity

F = flow rate


r = capillary radius

l = length

Measured Pressures in 50 μm ID Tubing




Measured Pressures in 25 μm ID Tubing




Influence of Mobile Phase on System Pressures

Backpressure in 25 μ m ID tubing using MeOH/water mixtures



A mixture of 60% MeOH in water has a much higher viscosity than either MeOH or water alone.

A Typical Split Flow System



ProteCol™ Low Dispersion Splitter

Control the sample flow to the column by changing the length of the bypass tube.


Can be used:

- Before the injector
- Between injector and column(s)
- Between column(s) and detector
- Extremely small internal dimensions and volume allow sample to pass through with minimal dispersion.




Splitter features a 2 μ m integral filter at the inlet.

[Splitter Sectioned Drawing




Splitter Section Close-up



Repeatability of Column and Splitter Performance


	Pyridine	MetBenz	Phenetole	Naphthalene
1st inj.	1.65	2.34	2.91	3.67
2nd inj.	1.66	2.34	2.92	3.68
3rd inj.	1.67	2.34	2.93	3.69
4th inj.	1.65	2.33	2.92	3.68
5th inj.	1.65	2.32	2.90	3.66
6th inj.	1.67	2.34	2.93	3.69
7th inj.	1.66	2.33	2.90	3.64
8th inj.	1.66	2.33	2.91	3.66
9th inj.	1.66	2.33	2.90	3.64
10th inj.	1.65	2.31	2.89	3.63
11th inj.	1.67	2.33	2.91	3.67
12th inj.	1.66	2.32	2.89	3.63
13th inj.	1.67	2.34	2.92	3.67
14th inj.	1.65	2.32	2.90	3.65
15th inj.	1.65	2.30	2.87	3.61
16th inj.	1.66	2.34	2.91	3.67
17th inj.	1.65	2.32	2.91	3.65
18th inj.	1.65	2.31	2.89	3.64
19th inj.	1.67	2.35	2.93	3.69
20th inj.	1.65	2.32	2.90	3.65
average	1.66	2.33	2.91	3.66
StdDev	0.01	0.01	0.02	0.02
%RSD	0.50	0.55	0.54	0.62




Column: 100mm x 300 μ m ProteCol C18, 3 μ m, 300Å
 Mobile Phase: 60%AcCN:40% water
 Flow Rate: 3.8 μ l/min
 Wavelength: 254 nm
 Flow splitting after injection valve



Overlaid Chromatograms from the Repeatability Experiment




Tryptic Digest of Cytochrome C




Robert Moritz, Joint Proteomics Laboratory
Ludwig Institute For Cancer Research & The Walter and Eliza Hall Institute of Medical Research
Parkville, Victoria, Australia 3050



Tryptic Digest of Ovalbumin



Sensitivity for Ovalbumin Digest




Comparison of Sensitivities

Ovalbumin Load	"Market Leader" Column		SGE ProteCol	
	No of Peptides	% Coverage	No of Peptides	% Coverage
1 fmol	X	X	6	26.0
8 fmol	5	20.8	8	30.1
32 fmol	6	30.9	10	41.3
125 fmol	8	35.1	13	48.6
500 fmol	9	35.8	15	51.2
1 pmol	12	42.9	17	51.2

Column: 75 μ m x 100mm, C18, 3 μ m, 300 \AA
Flow Rate: 150nL/min

Column: 150 μ m x 100mm, C18, 3 μ m, 300 \AA
Flow Rate: 550nL/min

ProteCol™ MicroFlow Meter



Specifications and features:

- Measures flow directly in a precisely calibrated syringe barrel
- Perfect companion to the flow splitter
- Two flow ranges available :
0.2 – 6 µl/min and 5 - 50 µl/min

[Conclusions

- Peak dispersion or bandspreading lowers HPLC efficiency, resolution and mass sensitivity.
- Extremely small peak volumes in capillary LC creates a special challenge for making column connections.
- Dispersion outside of the column bed must be strictly minimized in order to separate complex mixtures with high resolution.
- Wide acceptance of capillary LC and LC/MS techniques will require continuing improvement in system, column and accessory designs.

References

- Small-Bore Liquid Chromatography Columns, R. P. W. Scott, ed., Wiley, NY (1984).
- R. P. W. Scott, Extra Column Dispersion, Chrom-Ed Book Series.
- S. R. Bakalyar, et. al., Rheodyne Technical Notes 9, 1988.
- R. Majors, LCGC North America, Volume 21 Number 12, 1124-1133, December 2003.
- SGE unpublished report on performance of PEEK and PEEKsil[®] tubing under various operating conditions.
- Poster, “Performance of Chemically-Fabricated Frits in Packed Capillary LC Columns”, Pittcon 2004, Chicago, IL.
- Poster, “Nano-Flow LC-MS Columns for the Analysis of Oligosaccharides”, Pittcon 2004, Chicago, IL.
- Poster, “Multidimensional Approach to Capillary LC”, Pittcon 2004, Chicago, IL.
- Proteome Systems, unpublished proteomics LC-MS data.